

Selection of solvent of the probe used by the Microdrop Dispenser Heads

solvent
acetophenone
anisole
benzylalcohol
butylacetate
cyclohexanol
diethylene glycol
1,4-dioxane
DMA
DMF
DMSO
n-dodecane
ethylbenzene
ethylene glycol
hexylbenzene
hexyliodide
methyl-benzoate
N-methyl-pyrrolidone
1-octanol
1,3-propanediol
propylene carbonate
toluene
water
p-xylene

04/2016-e

Please contact us if you miss a solvent.



Advancing the Art of **Microdispensing**

microdrop Technologies GmbH

Tycho-Brahe-Kehre 1
22844 Norderstedt / Germany

Phone: +49 (0)40 53 53 83-0
Fax: +49 (0)40 53 53 83-24

info@microdrop.de
www.microdrop.de